

Gel Electrophoresis buffers and solutions

50x TAE Buffer – 1 L

Add in following order

683 ml DI H₂O
242.0 g Tris base
18.61 g Na₂EDTA•2H₂O
57.1 mL acetic acid
mix

To make 1x TAE – 1L
20 mL 50x TAE
980 mL DI H₂O

10x TBE buffer

Add to 1 L bottle

108 g Tris Base
55 g Boric acid
40 mL of 0.5 M EDTA (pH 8)

Fill to 1L with DI H₂O

5x Tris-Glycine-SDS buffer

Add to 1 L bottle

900 mL DI H₂O
15.1 g Tris Base
94 g glycine
50 mL 10% SDS

Fill to 1 L DI H₂O

6X loading dye (TAE)

add to 1.5 ml eppendorf tube

600 µL 50% v/v glycerol
120 µL 50x TAE
180 µL DI H₂O
100 µL Promega 6X loading dye

6X loading dye (TBE)

add to 1.5 ml eppendorf tube

600 µL 50% v/v glycerol
300 µL 10x TBE
100 µL Promega 6X loading dye

DNA Markers

Dilute marker to 50 µg/mL

add 20 % volume of 6X loading dye

Load 12 µL (0.5 ug) per lane

Ethidium Bromide (EtBr) 10mg/ml

add to 1.5 eppendorf tube

10 mg EtBr
990 µL DI H₂O