

## Diode pumped infrared 1064 nm CrystaLaser®

High reliability High stability Single longitudinal mode Ultra-compact


### SPECIFICATIONS

Wavelength	1064 nm	
Laser crystal	Nd:YAG/Nd:YVO <sub>4</sub>	
Laser version	Standard version	SLM version
Model No. and output power	CL1064-3W0 3W CL1064-2W5 2.5 W CL1064-2W0 2 W CL1064-1W5 1.5 W CL1064-1W0 1 W CL1064-800 800 mW CL1064-500 500 mW CL1064-300 300 mW CL1064-100 100 mW	CL1064-1W0-S 1 W CL1064-700-S 700 mW CL1064-500-S 500 mW CL1064-400-S 400 mW CL1064-300-S 300 mW CL1064-100-S 100 mW CL1064-050-S 50 mW CL1064-010-S 10 mW
Longitudinal mode spectrum	Multiple	Single; frequency tuneable option available
Coherence length	1 - 5 mm	> 300 meters
Linewidth	0.3 nm, nominal	< 10 <sup>-5</sup> nm
Modulation options	TTL 1kHz, 3 kHz, 10 kHz, 50kHz Analog control, Analog modulation	TTL on/off option
Output stability (rms)	< 1.5% over 24 hours; ultra-stable option: 0.5% or 0.25% over 24 hours	
Transverse mode	Single mode, TEM <sub>00</sub> , M <sup>2</sup> = 1.0 - 1.2	
Beam diameter (1/e <sup>2</sup> )	0.4 mm	Option: 2X - 5X laser beam expanders are available to change beam size and divergence e.g. 0.8 mm beam with 2 mrad divergence; or 2 mm beam with 0.8 mrad divergence
Beam divergence	4 mrad	
Beam pointing stability	< 0.005 mrad/°C	
Polarization ratio	100:1; Horizontal linear polarization	
Operation temperature	10 °C to 35 °C	
Warm-up time	< 1 minute	
Operation voltage	90-250 VAC or 12VDC OEM version 11-16VDC; Typical power consumption: 10 W	
Dimensions and weight of laser head	3x3x12 cm <sup>3</sup> (1.2x1.2x4.7 in <sup>3</sup> ) with a fixed 6 mm thick plate, 0.3kg (0.6 lb) 7x3.6x18.5 cm <sup>3</sup> (2.8x1.4x7.3 in <sup>3</sup> ), 0.6 kg (1.3 lb) for lasers ≥ 2 W, type 2 size	
Dimensions and weight of power supply	AC: 5x14x15 cm <sup>3</sup> (2x5.5x6 in <sup>3</sup> ), 0.5 kg (1 lb) DC: 8.5x3.5x12.7 cm <sup>3</sup> (3.3x1.4x5 in <sup>3</sup> ), 0.2 kg (0.5 lb); option CL-2005 power display AC: 5 x 16 x 15 cm <sup>3</sup>	



**Fiber coupling option:** Single mode and multi-mode fiber delivery.

**CL-2005 power supply option:** Power adjustable and power digital display.

**Made in the USA**




### Laser Head All dimensions in mm



Mechanical dimensions of the laser head and power supply of infrared laser are the same as the green laser.



Case I: In case, beam waist  $\omega_0$  is 5cm from the output window; outside of the laser; Our cavity should look like this.



We know beam waist at the output window  $M_2$   $\omega(z=-5\text{cm}) = 0.2\text{mm}$

Divergence  $\theta = 4 \times 10^{-3} \text{ rad}$

$$\lambda_0 = 1064 \times 10^{-6} \text{ mm}$$

The Beam waist at  $z=0$  is:  $\theta = \frac{2\lambda}{\pi \omega_0} \Rightarrow \omega_0 = \frac{2\lambda}{\pi \theta}$

$$\omega_0 = 0.169 \text{ mm} \approx \boxed{\omega_0 = 0.17 \text{ mm}}$$


To check, if this agrees with the specified beam waist at  $z=-50\text{mm}$

$$\omega(z=-50\text{mm}) = \omega_0 \sqrt{1 + \left(\frac{z}{z_0}\right)^2}$$

$$\text{where } z_0 = \frac{\pi \omega_0^2}{\lambda_0} \Rightarrow \boxed{z_0 = 85 \text{ mm}}$$

$$\boxed{\omega(z=-50\text{mm}) = 0.197 \text{ mm} \approx 0.2 \text{ mm}} \quad \begin{array}{c} \downarrow \downarrow \\ \text{It AgRees} \end{array}$$

Now with 4x option:



magnification  $M = 4$

waist  $\omega_0' = M\omega_0$

Divergence  $\theta_0' = \frac{\theta_0}{M}$   $\theta_0 = 4 \text{ mrad}$

new  $Z_0' = M^2(2Z_0)$

$$Z_0' = 1360 \text{ mm}$$

Beam waist at the ~~output~~ output window  $\omega_1' = \omega_1 \times M$   $\omega_1 = 0.2 \text{ mm}$

new beam waist at new output window  $\omega_1' = 0.2 \times 4 = 0.8 \text{ mm}$

Divergence  $\theta_1' = \frac{4}{4} = 1 \text{ mrad}$


Now the new beam waist at  $Z=0$   $\omega_0' = \frac{2\lambda}{\pi\theta_1'} = 0.68 \text{ mm}$

$$\omega_0' = 0.68 \text{ mm}$$

also new  $Z_0' = \frac{\pi\omega_0'^2}{\lambda_0} = 1360 \text{ mm}$

Now the distance  $d'$  where the new beam waist exists from the output window

$$\omega_1' = \omega_0' \sqrt{1 + \left(\frac{Z}{Z_0'}\right)^2} \quad [Z = d']$$




$$\frac{\omega'_0}{\omega'_1} = \frac{1}{\left[1 + \left(\frac{\pi \omega'_1{}^2}{\lambda R'}\right)^2\right]^{1/2}}$$

$$\Rightarrow \boxed{R' = 3047 \text{ mm}}$$

$$\text{now } d' = \frac{R'}{1 + \left(\frac{\lambda R'}{\pi \omega'_1{}^2}\right)^2} \Rightarrow \boxed{d' = 846 \text{ mm}}$$

Case II In case, the beam waist is 5cm from the output window; Inside of the laser cavity! Our cavity should look like this



~~Every~~ Beam waist,  $z_0$  & divergence will remain <sup>the</sup> same  $w_0 = 17\text{mm}$   
 $z_0 = 85\text{mm}$   
 $\theta = 4\text{mrad}$

Curvature of the output window  $R_2 = 195\text{mm}$