



System Engineering

20.109(F10)
M2D5 lecture
10.28.10


New tools for reliable engineering of complex biological systems




Personal Genomics




Synthetic Genomics



Beyond the C-dog



Genetic screen









K+ Library Variations (in red)

Cph1/EnvZ	A553	G554	V555	S556	H557
EnvZ	A239T	G240E	V241G	S242D	H243A
wt seq	GCG	GGG	GTA	AGT	CAC
oligo seq	RNS	RNS	RNS	RNS	SNW
poss aa	Val	Val	Val	Val	Val
	Ala	Ala	Ala	Ala	Ala
	Asp	Asp	Asp	Asp	Asp
	Glu	Glu	Glu	Glu	Glu
	Gly	Gly	Gly	Gly	Gly
	Ile	Ile	Ile	Ile	Leu
	Met	Met	Met	Met	Pro
	Thr	Thr	Thr	Thr	His
	Asn	Asn	Asn	Asn	Gln
	Lys	Lys	Lys	Lys	Arg
	Ser	Ser	Ser	Ser	
	Arg	Arg	Arg		

NOTE: no stop codons should be in mix

R = G, A
 N = G, A, T, C
 S = G, C
 W = T, A





20.109(F10): Laboratory Fundamentals of Biological Engineering

(1) Sequencing K+ candidates
Miniprep DNA as you did in Module 1
(Soln I, Soln II, Soln III, EtOH, wash, dry)

Resuspend pellets in 40 ul of water

For sequencing, mix:

- 2 ul plasmid DNA
- 6.4 ul of a 1:100 dilution of the primer NO289
- 15.6 ul sterile water

20.109(F10): Laboratory Fundamentals of Biological Engineering

(2) Also check your DNA by digest

- Provide 10 ul of each
- Teaching faculty will cut with

Nde I
Xba I

to run on agarose gel and post

