


# Microrheology of Complex Fluid



❑ Rheology: Science of the deformation & flow of matter

❑ Microrheology


- Microscopic scale samples
- Micrometer lengths



Complex shear modulus  $G^*(\omega)$


$$\sigma = G^* \varepsilon$$

- $G^*(\omega) = G'(\omega) + j G''(\omega)$
- Solid vs. fluid
- Resistance to deformation



ciks.cbt.nist.gov

Storage modulus  $G'$   
Energy storage  
Elasticity ~ Solid



ma.man.ac.uk

Loss modulus  $G''$   
Energy dissipation  
Viscosity ~ Fluid

# High Frequency Microrheology Measurement


## Active Method:

Magnetic microrheometer – Baush, BJ 1998  
Huang, BJ 2002


## Passive Method:

Single particle tracking – Mason, PRL 1995  
Yamada, BJ 2000  
Multiple particle tracking – Crocker, PRL 2000

# Magnetic Microrheology



**Magnetic Microrheology**



**5 sec Step Response**

# Basic Physics of Magnetic Microrheometer



Ferromagnetic particle

$$\mathbf{F} = \frac{1}{2} m_0 \nabla (\mathbf{m} \cdot \mathbf{H})$$

Particles cluster together!  
Doesn't work!

Paramagnetic particle – no permanent magnetic moment

$$\mathbf{F} = m_0 c V \nabla (\mathbf{H} \cdot \mathbf{H})$$


$\chi$  is susceptibility

V is volume

Note: (1) force depends on volume of particle  
(5 micron bead provide 125x more force)  
(2) force depends on magnetic field GRADIENT


# Magnetic manipulation in 3D

Top View



ST: stage  
PH: post holder  
CO: coil  
(400 turn/cm)  
PO: pole  
SA: sample chamber  
OB: objectives  
(100x, 1.0 n.a.  
water;  
20x, 0.5 n.a  
all reflecting)

Side View



ST


\*Lower Force  
nN level

\*3D

\*Uniform gradient

Amblad, RSI 1996  
Huang, BJ 2002

# Magnetic manipulation in 1D



\*High force  
>10 nN


\*Field non-uniform  
Needs careful  
alignment of tip  
to within microns

\*1D

Baush, BJ 1998

The bandwidth of ALL magnetic microrheometer is limited by the inductance of the eletromagnet to about kiloHertz

# Magnetic Rheometer Requires Calibration





Baush, BJ 1998

# Mag Rheometer Experimental Results



Baush, BJ  
1998




Transient responses allow fitting to micro-mechanical model

Problem – Magnetic bead rolling

Solution – Injection, Endocytosis Modeling (Karcher BJ 2003)

# Model Strain Field Distribution



Baush, BJ 1998