


Advanced Fluorescence Microscopy




Spectral Resolved Microscopy

Lifetime Resolved Microscopy


Foster Resonance Energy Transfer (FRET) Microscopy

Fluorescence



Jablonski Diagram




Spectral Resolved Fluorescence Microscopy




In collaboration with Dr. C. Buehler, PSI






Spectrally Resolved 2-P 3-D Imaging




3 Color Latex Beads
(10 Hz frame rate)



Ex vivo human dermis
(0.2 Hz frame rate)

Fluorescence

Radiative & Non-Radiative Decay I



Fluorescence

Radiative & Non-Radiative Decay II

Fluorescence lifetime

$$t = \frac{1}{\Gamma + k}$$

Fluorescence Natural lifetime

$$t_0 = \frac{1}{\Gamma}$$

Fluorescence

Intensity Decay

Single Exponential Decay

$$\frac{dN_e}{dt} = -(k + \Gamma)N_e$$

$$F = F_0 e^{-(k+\Gamma)t} = F_0 e^{-t/\tau}$$

Fluorescence

Intensity Decay II

Multiple Exponential Decay

$$\frac{dN_e}{dt} = -\left(\frac{1}{t_1} + \frac{1}{t_2} + \dots\right)N_e$$

$$F = F_0 e^{-\Sigma t/t_i}$$

Fluorescence

Lifetime of Fluorophores

POPOP: 1.3 ns

Fluorescein: 4 ns

GFP: 2.5 ns

CdSe/ZnS Q-dot: 25-50 ns




Terbium, Rubidium compounds: >500ns

EtBr (free): 1.8 ns


EtBr (bound): 23 ns

Time domain

Time Correlated Single Photon Counting




Time-Correlated Single Photon Counting



Lifetime imaging and biological functions


- (1) Distinguish cellular organelles by multiple lifetime imaging

Lifetime Sections




- (2) Monitor metabolite concentration (Ca, pH etc)

Intensity




Lifetime




- (3) Monitor cellular processes such as proteolytic activity

Intensity




Lifetime




Lifetime Image of ECM in Dermis

Mean lifetime Image




Lifetime decay at a selected point



Zero lifetime component:
Collagen 2nd harmonic
emission

ns lifetime component:
Elastin fluorescence

Resolving Collagen & Elastin Fraction in the Dermis





Fraction of SHG (Collagen)


Foster (Fluorescence) Resonance Energy Transfer

Focal Adhesion Complex – Paxillin-Focal Adhesion Kinase Interaction


- *Double negative mutants of FAK is fatal. Tumor cells over expresses FAK
- *FAK controls cell spreading, migration, and survival
- *Focal adhesion localization results in auto-phosphorylation of FAK
- *FAK phosphorylation is linked to cytoskeletal remodeling (Src, Rac 1)




(1) Paxillin-FAT in endothelial cells




GFP-Paxillin




FAT-mCherry



Spectral overlap



Foster (Fluorescence) Resonance Energy Transfer



Dipole - dipole interaction

r^6 dependence

Efficiency

50% energy transfer

Förster distance

$R_0 = 40$ to 70 \AA

Decrease excited-state lifetime

Fixed geometry yields constant


FRET parameters

$$E = \frac{R_0^6}{R_0^6 + r^6} = 1 - \frac{t_{DA}}{t_D}$$


$$R_0^6 = \frac{9000 \ln(10) k^2 f_D}{128 p^5 N_A n^4} J$$

$$\text{where, } J = \frac{\int F_D(I) e_A(I) I^4 dI}{\int F_D(I)}$$

FRET & Spectral Bleedthrough









Elangovan, Methods, 2002



FRET Image: Are we ok now if we see acceptor fluorescence?




Spectral Resolved FRET




Layout for image acquisition (one- or two-photon excitation) for precision FRET (PFRET) data analysis^a

Wavelength excitation	Fluorophore	Emission donor images	Emission acceptor images
λ_D	D	a	b
	A		c
	D + A	e	f
λ_A	A		d
	D + A		g

(1) Lifetime Resolved FRET





“Quantify” Signaling Pathway Using t-FRET



$$I_i^{\text{model}}(t) = \int_0^t G(t-T) \times c_{2i+1} \left(c_{2i+2} \exp\left(-\frac{T}{c_1}\right) + (1 - c_{2i+2}) \exp\left(-\frac{T}{c_2}\right) \right) dT.$$

(1) Resolving bound paxillin fractions In PAC and cytosol



FRET ratio

0.17


0.06

0.003

0

0

Increasing intensity



FRET ratio

0.54 ← cytoplasm

0.84 ← majority FA

1

0.65