

Measurement of higher brain function in humans

Michael S. Beauchamp, Ph.D.



Assistant Professor
Department of
Neurobiology and
Anatomy



Adjunct Assistant
Professor
Departments of
Psychology and
Bioengineering



Adjunct Assistant
Professor
Department of
Neuroscience

Michael.S.Beauchamp@uth.tmc.edu

*Higher Brain Functions, Mariella De Biasi, Course
Director*

April 26/28, 2011


Introduction to fMRI

Cross-listed as:

Rice Psychology 671

BCM GS-NE-439

Textbook



Functional Magnetic Resonance Imaging by
Scott A. Huettel, Allen W. Song, and Gregory
McCarthy (Hardcover - April 1, 2004)

Why is neuroimaging difficult?



Palpation/Sensation

- skull: can't feel much from outside (sorry, ultrasound)
- no nerves in brain: can't feel much from inside

Alphabet Soup

CT

MRI

SPECT


EEG

fMRI


PET

MEG


Taxonomy of Neuroimaging Methods



Taxonomy of Neuroimaging Methods




Taxonomy of Neuroimaging Methods



X-Ray (radiograph)—bones, not brain!
(tissue density)



Taxonomy of Neuroimaging Methods




Arteriogram (a.k.a. Angiogram)


Basic principle: Inject contrast agent (dye) that is radio-opaque i.e. iodine containing agents



Taxonomy of Neuroimaging Methods



CT (computed tomography)



CT (computed tomography)

Pros:

Widely available

Very fast to collect whole-head images

(one slice in < 1 ms; whole head in ~ seconds; whole exam in 10 minutes)

Somewhat cheaper than MRI
(~\$500 vs ~\$1000 for MRI)

Less hassle (few contraindications)

Best for an emergency


Cons:

Exposure to ionizing radiation

(increased risk of cancer)


can require a contrast agent—for brain, injectable iodine compound

Poor tissue contrast
not versatile




CT (computed tomography)

Basic Principle: rotate machinery to take multiple x-rays with different paths through the body



Common clinical use: stroke

Patient presents with stroke



Give tPA, dissolve clot,
Blood flow restored,
Patient recovers

Give tPA, prevent clotting,
Patient dies of massive bleed

Common clinical use: stroke

Patient presents with stroke



80%

ischemic



20%

hemorrhagic



CT (computed tomography)

Pros:

Widely available

Very fast to collect whole-head images

(one slice in < 1 ms; whole head in ~ seconds; whole exam in 10 minutes)

Somewhat cheaper than MRI
(~\$500 vs ~\$1000 for MRI)

Less hassle (few contraindications)

Best for an emergency


Cons:


Exposure to ionizing radiation

(increased risk of cancer)


can require a contrast agent—for brain,
injectable iodine compound

Poor tissue contrast
not versatile






“Nuclear Medicine”



“Nuclear Medicine”: PET/SPECT



Basic principle:

Inject radioactive isotope attached to metabolic compound (Oxygen, Glucose, etc.)

Wait for it to decay, look for radioactive decay products

PET/SPECT

Pros:

Fairly cheap (~ \$1000)

Shows function
(metabolism)

Cons:

Exposure to ionizing
radiation (increased risk
of cancer)

~1 year of background
radiation

Low-resolution


Slow to very slow
not versatile

SPECT

Single Photon Emission Computed Tomography




“Nuclear Medicine”:
~unclear medicine



SPECT

Single Photon Emission Computed Tomography



Radionuclides

Single-Gamma emitting

^{99m}Tc , ^{123}I , ^{67}Ga , ^{111}In

PET: Positron Emission Tomography




Basic principle:

Inject radioactive isotope attached to important metabolic compound (Oxygen, Glucose)

Wait for it to decay.

Pick up **two** particles going in opposite directions—improves spatial resolution



- PET uses beta-plus-emitting radionuclides such as C-11, N-13, O-15, and F-18 which annihilate into two 511-keV photons that travel in opposite directions.


Developments in PET

■ Development of new radiotracers

[11C]DTBZ

VMAT2

~dopamine, serotonin




DVB

[11C] Flumazenil

(benzodiazepine receptor antagonist


~ GABA-A)




DVR

[11C]PMP


Substrate for AChE



2 –photons instead of one → better resolution



Sample Clinical Application: Alzheimer's Diagnosis




PiB PET SCANS

- Compounds that bind to AD plaques
- Conditional approval by FDA (Jan 2011)
- C. M. Clark et al. J. Am. Med. Assoc. **305**, 275–283; 2011



University of Pittsburgh
PET Amyloid Imaging Group



MRI: Magnetic Resonance Imaging

Pros:

incredible images

advancing very rapidly

extremely high resolution

extremely versatile

NO ionizing radiation

Cons:

moderately expensive (~
\$1000)

complex



some contraindications

Can require injection of
contrast agents (Gadolinium/
Iron compounds)

Imaging Techniques: MRI

- the MR scanner is a giant magnet: 1.5T, 3T, 7T, 9T






Contraindications I

- Ferrous metal in body




Contraindications II

- Cochlear implants (always)
- Maybe: pacemakers, vagal nerve stimulators, old (> 20 year) surgical implants



Clinical Applications



BCM HNL



- 5 research-dedicated 3 tesla scanners made by Siemens
- 2 head-only, 3 full-body




MRI (Magnetic Resonance Imaging)

Basic principle: uses radio waves to interrogate protons in water molecules in the brain



MRI (Magnetic Resonance Imaging)


Basic principle: uses radio waves to interrogate protons in water molecules in the brain



128 MHz at 3T (~ FM Radio)



We listen to these radio waves with an “RF coil” (radio antenna)



H₂O



H₂O




H₂O




B_0 = Giant Field Produced by Giant Magnet

- Purpose is to align H protons in H_2O (little magnets)




[Main magnet and some of its lines of force]

[Little magnets lining up with external lines of force]



◆ Small B_0 produces small net magnetization M

◆ Thermal motions try to randomize alignment of proton magnets





◆ Larger B_0 produces larger net magnetization M , lined up with B_0

◆ Reality check:
0.0003% of protons aligned per Tesla of B_0

Precession of Magnetization M



- Magnetic field causes M to rotate (or *precess*) about the direction of B at a frequency proportional to the size of B — 42 million times per second (42 MHz), per Tesla of B



- ◆ N.B.: part of M parallel to B (M_z) does not precess


B_1 = Excitation (Transmitted) Radio Frequency (RF) Field

- Left alone, M will align itself with B in about 2–3 s
- So don't leave it alone: apply (transmit) a magnetic field B_1 that fluctuates at the precession frequency and points perpendicular to B_0




Relaxation: Nothing Lasts Forever

- In absence of external B_1 , M will go back to being aligned with static field B_0 — this is called *relaxation*
- T2: Part of M perpendicular to B_0 shrinks [M_{xy}]
 - This part of M is called *transverse magnetization*
 - It provides the detectable RF signal
- T1: Part of M parallel to B_0 grows back [M_z]
 - This part of M is called *longitudinal magnetization*
 - Not directly detectable, but is converted into transverse magnetization by externally applied B_1




Ants on a Pole analogy—two different
physical properties


T1 recovery




T2 decay



T1WI




T2WI




Basics of MR

- T1 ~ Longitudinal
Magnetization/Relaxation
“hi-resolution, normal anatomy”




T2 ~ Transverse Magnetization/
Relaxation
“pathology”—water content




Basics of MR—contrast agents

- Gadolinium
- can be injected to enhance contrast (usually in T1 images); hastens T1 recovery making image brighter



Clinical Applications: Multiple Sclerosis

- T2 is best for seeing white matter abnormalities




Bonus—Shelf Exam Question

A 53-year-old woman dies 4 days after an automobile collision. She sustained multiple injuries including a femoral fracture. Widespread petechiae are found in the cerebral white matter at autopsy. Which of the following is the most likely cause of these findings?

- (A) Acute respiratory distress syndrome
- (B) Contrecoup injury
- (C) Fat embolization
- (D) Septicemia
- (E) Subdural hematoma

Fat emboli (FLAIR T2)



Scanner as computer...


- pulse sequences are software

CT vs MRI

TV with one channel (not that much good on)


vs.

Computer that can run lots of programs (takes longer to boot up, but get more for it)



MRI vs. fMRI

MRI studies brain anatomy.



Functional MRI (fMRI)
studies brain function.



An “anonymous” grad student





MRI vs. fMRI

high resolution
(1 mm)


MRI



one image

fMRI

low resolution
(~3 mm but can be better)



many images
(e.g., every 2 sec for 5 mins)


fMRI

Blood Oxygenation Level Dependent (BOLD) signal
indirect measure of neural activity

↑ neural activity → ↑ blood oxygen → ↑ fMRI signal



An “anonymous” M3





Metabolism

- Brain uses ~20% of total body oxygen
(even though only 1-2% of total body mass)



Complex mechanism for regulating cerebral
blood flow to ensure adequate oxygen supply

Deoxygenated blood attenuates T_2^* -weighted MR images



MRI volume element


decrease of venous dHb during increased perfusion:




Human 3T used for BOLD fMRI




One Fast Image



Overlay on Anatomy






Graphs vs. time of 3×3 voxel region




X: 29
Y: 29 Grid: 20 Scale: 5 pix/datum
Z: 6 Num: 68 Base: separate


68 points in time 5 s apart; 16 slices of 64×64 images








Diffusion Weighted Imaging—earliest post-stroke diagnosis





MR Angiography




Magnetic Resonance Spectroscopy


B




reduced GABA in panic disorder

Cortical Surface Models—AD/PD








(a)




Excitatory synapse



Inhibitory synapse



(b)



EEG/MEG

Pros:

Completely non-invasive—**only modality that NEVER requires injection of a contrast agent**

no ionizing radiation

direct measurement of neuronal activity


Cons:

not really neuroimaging

Limited clinical utility:


EEG—sleep studies

MEG--epilepsy




EEG

basic principle: electrodes on scalp surface record summed electrical activity (mainly synaptic) of many neurons




MEG Scanner

basic principle: sensors near scalp surface record summed magnetic field resulting from electrical activity (mainly synaptic) of many neurons




EEG/MEG Results

EEG Activation Map



MEG Activation Map



Taxonomy of Neuroimaging Methods

