

- Announcements

- No lab next week

- OH: NK by apt. ANS by apt.

- Lab Quiz

- Pre-lab Lecture



- ❖ ECD Preparation

- ❖ Intro to TEM

- ❖ Today in Lab

electro deposition


Preparing ECD, part 1



* don't dip tweezers into liquid

TEM

- Very high resolution – why? e^- has much lower λ than light
- EM lens to focus
- Sample preparation
 - very thin, under vacuum
 - can't image *in situ* bio.
- cryoTEM *frozen, more in situ-like*
- S(canning)TEM *-today*
 - scanning allows mapping
 - = \uparrow efficiency, \uparrow contrast
- Unscattered e^- visualized according to sample density *"brightfield"*



Today in Lab

- Three groups to TEM first
 - yellow, blue, purple
 - use *your* phage, with equal volume IrOx
 - TEM with Mark Allen, room 13-1012
- Two groups set up electrodeposition first
 - green, pink
 - use *stock* phage, 1:50 with IrOx
↳ high PFU
- Then switch off