

UNIVERSITY OF ILLINOIS
AT URBANA-CHAMPAIGN

Illinois iGEM 2009

Bacterial Decoder



illinois.edu



Steve
Waltersdorf

Sonam Tanna

David Reif

Davel Patel

Dave Korenchan

Harsh Shah

Lucy
McCommas

Graham
Helmberg

Heather
August

Francis
Lee

Meagan
Musselman

Phil
Bell


Our Goal

“To create a universal, **modular**, and dynamic bacterial **decoder** to be used as a tool in a variety of **cellular logic** applications, while also contributing novel **BioBricks** to the Registry of Standard Biological Parts.”




The Bacterial Decoder

- A logic system that identifies combinations of specific inputs to express outputs
 - 2:4 3:8 4:16 $n:2^n$
- Involves AND gates, NOR gates, and inverters
 - Made from regulation on translational and transcriptional levels



IPTG	aTc	FP
0	0	C
0	1	R
1	0	Y
1	1	G




A Genetic "if" Statement

```
if (IPTG && aTc) {  
    return GFP;  
}
```

```
else if (IPTG && !aTc) {  
    return YFP;  
}
```

```
else if (!IPTG && aTc) {  
    return RFP;  
}
```

```
else if (!IPTG && !aTc) {  
    return CFP;  
}
```




Use and Future Applications


Sensing for
multiple water
contaminants



Complex logic
for cellular
computer



Medical
diagnostics and
treatment




Approach

Gene regulation using transcription factors and small RNAs.





Small RNAs (sRNAs)

- Can repress protein synthesis
- Mechanisms of *trans*-encoded Hfq binding inhibitory sRNAs:
 - RBS occlusion
 - Degradation



sRNA Advantages




- Quick response times are ideal for state changes
- sRNA regulatory mechanisms are cost effective



Implementing sRNA in a decoder

- Small RNAs are used in concert with transcription factors to regulate gene expression.
- Created a NOR gate with TF and Small RNAs.



Decoder Schematic








Diagram of Method






Transform plasmids
into cells

E.Coli Top
10 cells on
LB/cat




E.Coli Top
10 F' cells
on LB/amp




Co-transform

E. Coli Top
10 F' Cells




Genetic Circuit



sRNA preparation

- sRNAs should have a blunt 5' end
 - Incompatible with the standard BB format
- Under standard lac promoter
 - BBa_R0010
- Typically expressed alone under a promoter




Target sequence preparation

- RBS containing target sequences are fused upstream of the GFP reporter
- Standard BioBrick “prefixes and suffixes” would not allow for an in frame fusion.

Standard BioBrick assembly


...ACtactagagCA...



8bp → Frame shift

Freiburg Assembly Standard

...ACtactagCA...




6bp → No frame shift



MicF regulation

MicF sRNA/Target regulation



Submitted Parts

sRNA gene BioBricks

MicC

MicF

MicA

SgrS

GcvB

DicF

5' Target Sequence BioBricks

OmpC

OmpF

OmpA

Ptsg

Oppa & Dppa

Ftsz

galK



Status and Conclusion

- Biology is hard.
- We created several of the constructs for the decoder, but didn't have time to get it in working order
- We are submitting 14 novel BioBricks - the first sRNAs!



Advisors



Dr. Ido Golding



Courtney Evans



Dr. Chris Rao



Acknowledgments



MAX-PLANCK-GESELLSCHAFT

The Vogel Lab at the Max
Planck Institute



ENGINEERING AT ILLINOIS



illinois.edu



QUESTIONS?

