

Detection of Mutations in EGFR in Circulating Lung-Cancer Cells

Shyamala Maheswaran, Lecia V. Sequist, Sunitha Nagrath, et al.

New England Journal of Medicine 2008

20.309 Fall 2010


Helen Chen, Jennifer Lai

Cancer Treatment: Drug Resistance

- Mutations acquired during treatment → drug resistance
 - Non-small-cell lung cancer: EGFR activating mutations
 - Tyrosine kinase inhibitors eg gefitinib and erlotinib
 - Relapse within 1 year: acquisition of secondary EGFR mutation T790M
- Invasive biopsies provide insufficient material for molec analysis and cannot be sampled repeatedly
- Alternative biomarkers in the blood
 - Cell-free DNA – insensitive detection of tumor growth
- Molecular characterization of circulating tumor cells (CTCs)


CTC-chip

- Isolate, quantify, analyze CTCs
- Blood flows past EpCAM-coated microposts
- Optimize laminar flow conditions → high yield and purity, reproducible
- Correlated with disease progression



CTC Identification and Enumeration by Fluorescence Microscopy

- Cells fixed on CTC-chip
- Fluorescence Staining
 - Hoeschst → DNA content
 - Cytokeratin → epithelial cells
 - CD45 → leukocytes
- Identify CTCs based on cell size, morphology, and fluorescence



Supplementary figure 1: CTC-chip analysis of blood specimens from NSCL patients

EGFR Mutational Analysis

- Allele-specific mutation detection by Scorpion Amplified Refractory Mutation System (SARMS) in RT-PCR
- Scorpion primer:
 - Quenched fluorophore, gene-specific probe region, and primer to WT or mutant allele
- Primers anneal to DNA
 - Fluorophore still quenched
- Denaturation
 - Self-association of Scorpion primer → dissociation of quencher




Image from Clinical Cancer Research July 2006;12: 3875


Results

- Identified circulating tumor cells
- Detected EGFR mutations in tumors
- Tumor genotyping over time

Identifying CTCs

- Patients:
 - 23 EGFR mutant tumors
 - 4 WT EGFR tumors
 - 10 mL blood sample
 - Radiographic measurements performed close to time of CTC analysis
- Median of 74 CTCs/mL (mean, 133; range, 5-771)
- Quantity of CTCs did not correlate with tumor volume as assessed by radiographic measurements

Detection of EGFR mutations in tumors







- Mutations:
 - Exon 19 mutations (“del”)
 - L858R missense mutation
 - T790M mutation -> confers resistance to gefitinib and erlotinib
- T790M mutation in pre-treatment tumor samples
 - 10/26 patients
 - High number of amplification rounds

- pre-treatment T790M mutation correlated with progression-free survival
 - 7.7 months vs. 16.5 months

CTC count over time

- Look at effect of treatment on CTCs






— CTCs per mL

— Tumor size

[] low allele frequency

Patient 9

B



Treatment-induced genotype shift

Significance

- CTC count may be a good measure of tumor metastasis and invasiveness
 - Better predictor than simply tumor size
- Useful for tumor sub-typing
 - Pre-existence of T790M mutation is a predictor of response to gefitinib/erlotinib treatment
- Track tumor genotypes over time
 - Cancer treatment is a selective process
 - Might aid medical treatment decisions

Limitations and Future Directions

- SARMS limited to known mutations
 - Cheaper/faster sequencing
- Not high-throughput
 - Combine with:
 - Microfluidic DNA isolation
 - Microfluidic PCR
- Application to other cancers with well-characterized mutations

Sources

- Maheswaran, *et al.*, "Detection of Mutations in EGFR in Circulating Lung-Cancer Cells" NEJM 2008.
- Whitcombe D, Theaker J, Guy SP, Brown T, Little S. Detection of PCR products using self-probing amplicons and fluorescence. Nat Biotechnol 1999;17:804-7.
- Nagrath S, Sequist LV, Maheswaran S, et al. Isolation of rare circulating tumour cells in cancer patients by microchip technology. Nature 2007;450:1235-9.

Validate SARMS assay

- Designed to detect several mutations in EGFR
 - Exon 19 mutations (“del”)
 - L858R missense mutation
 - T790M mutation -> confers resistance to gefitinib and erlotinib
- Compared SARMS and sequencing results
 - Identified the same mutation in 25/26 EGFR-mutant tumors




Image from Whitcombe, et al.1999


Comparison of CTC analysis and free plasma DNA analysis

- 18 patients
 - CTCs
 - Free plasma DNA
- 12 patients for sensitivity calculation
 - Primary tumor sample
 - CTCs
 - Free plasma DNA

| | EGFR mut | Sensitivity |
|--------|----------|-------------|
| CTCs | 17 | 92% |
| plasma | 7 | 33% |

SARMS

- Allele-specific primers compose of specific fluorophore, stem, gene-specific probe region, quencher, blocker molecule, and primer to WT or mutant
- PCR with Scorpion primers specific to WT or mutant alleles
- Allele-specific extension



Statistical Analysis

- Statistical analysis on
 - CTC quantity vs tumor burden
 - Mutations identified in different populations
 - Baseline T790M vs progression-free survival

Patients and Clinical Specimens

Table 1. Detection of Circulating Tumor Cells in Patients with Non–Small-Cell Lung Cancer.*

| Patient No. and <i>EGFR</i> Mutation Status | Sex | Age yr | Histologic Features | Time since Diagnosis mo | Previous Systemic Therapy† | Tumor Burden‡ | Circulating Tumor Cells§ |
|--|-----|-----------|---------------------|-------------------------------|-------------------------------|------------------|-----------------------------|
| | | | | | | cm | no. per ml |
| <i>EGFR</i> mutation present | | | | | | | |
| 1 | M | 58 | Adeno | 3.2 | None | 19.8 | 156 |
| 2 | M | 55 | Adeno | 14.4 | C, E | 2.4 | 50 |
| 3 | F | 66 | Adeno | 18.2 | G, C | 19.5 | 9 |
| 4 | M | 59 | Adeno/BAC | 20.7 | G | 2.0 | 771 |
| 5 | M | 57 | Adeno | 10.8 | C | 1.5 | 152 |
| 6 | F | 74 | Adeno | 18.3 | E | 9.8 | 5 |
| 7 | M | 64 | NSCLC | 13.1 | G, E, C | 4.6 | 196 |
| 8 | F | 70 | Adeno | 1.4 | None | 5.8 | 175 |
| 9 | F | 63 | Adeno | 0.9 | None | 28.6 | 143 |
| 10 | F | 66 | Adeno/BAC | 1.3 | None | 8.2 | 112 |
| 11 | F | 74 | Adeno/BAC | 4.8 | G | 4.5 | 74 |
| 12 | F | 62 | Adeno/BAC | 56.8 | G, E, O | 7.2 | 9 |
| 13 | M | 27 | Adeno | 9.9 | G | 4.1 | 47 |
| 14 | F | 55 | Adeno | 11 | G | 5.2 | 241 |
| 15 | F | 70 | Adeno/BAC | 54.7 | G, C | 30.2 | 31 |
| 16 | M | 53 | Adeno/BAC | 9.2 | E | 7.0 | 49 |
| 17 | F | 60 | Adeno/BAC | 97.0 | G | 4.3 | 103 |