

# *Our Field of Dreams*

## **Team Emblazon**

Diwash Acharya

Sushma Kalle

Rajiv Manjunath

Angie Nelson

Binu Maria Sacharias

Melissa Sommer


Tiffany Walters





# Original Plan

- Isolate DNA and RNA
  - Amplify using primers designed from mRNA sequence on GENBANK and that would add BioBrick prefix and suffix
- Sequence DNA and analyze for any introns
- Isolate vector and promoter from IGEM 2007 parts registry
- Insert new part and grow on naringenin enhanced media to produce blue colonies of *E. coli*




# DNA Isolation

- Used Hops DNA Isolation protocol from Dr. Schwekendiek's lab




# RNA Isolation

- Used Rneasy Kit from QIAGEN



# PCR Amplification


- Worked once, then failed repeatedly until it worked again



- Varied Amount of DNA (1  $\mu$ l, 2  $\mu$ l, 4  $\mu$ l, 8  $\mu$ l, 10  $\mu$ l) and temperature (53.7°C and 58°C)


# Gene Sequence

- >B01\_p450-F\_F1BBa\_528838.Seq Sequence #0 of 7 downloaded on Fri Nov 21 11:40:23 CST 2008
  - NNNNNNNNNNNNNNTNNCNNCTCTTCCTCTCTCTTTNCTCCCCTCCTACTCCGCCGCTTCCGCCGCCGCCGCCGCATCCTCCCCCTCCCTCCCGACCCCTCAACTTCCCCATCGTCGGCGCACTCCCTTCATCGGCCCATGCCCACTCAGCCTCGCCCTCCTCTCCCGCCGCTACGGCCCCATCATGTTCTCAAAATGGGCATCCGCCAAGTGGTCGTCGCCTCCTCCTCCTCCGCCGCACGCTCCTTTCTCAAACTCAGACTCCCGCTTCTCCGACCGCCCTCTCGACATCATCTCCAATCAAGTCAGCTACAACGGCCAGAACAGGGTCTTCGCCGACTACGGTCCCAAGTGAAGCNCCTCCGCAAA GTCTCCAATCTCCATCTCTTTGGCCCCAAGGCCATGTCCCGTTGGGCCGATGTGCGGCGGGAN GAAGCTTTCTCCATGTCCCACTTCCTGAAGAAACAGAGCGATTCCAAAAACCCTGTTTTGCTTTC CAACTTGCTGGTTTGTTCCATGGCGAATGTGATTGGGAGAATTGCAATGAGCAAGAGGGTGTTG GATGANGAGGGGAAGGAGGCNAACGAGTTTAAGGAGANANTNAANGANCTTTTGGNCGGGCA GNNNGCTTCNAATATTGNNGANTTGGTGCCNGCGANNANNNGGCTGGANCCTCNNGNNGNNA NGANNANNNGCCNGNNANTGAATCTNANGTTTGNNANGATGATANNTTAGTTCCTGGCGGAT CACGGTGACANTANAGGGNANCTGGCANGTGAATCCNGATCTGCTTGACCTTATTGTGGNTGA NANGANCCCCGGNNANAACGTNNANNNNCTTTTTCTNNANAACACCTCTANGTCTTCTTATCTTG TAATCAAANCTNNAATTCTNGATTTTANTTTNAATAAACGACTAATNNTNNTAATANNNNACTTTTT CTTNAANGANTCNTATNANNCTACATCTCNTCCTCCANAANTNNANNTANN
- >C01\_p450-R\_R1BBa\_528839.Seq Sequence #1 of 7 downloaded on Fri Nov 21 11:40:23 CST 2008
  - NNNNANNNNNNNNTTTTNNNNNNNNNGCCGCCGGCAGCCGCGGCCTCGCCGTCNCCAAGAGCGG CACCGCCTTCGGCAATACCAATCCCGGCCCTCCTCCATGTCCAGCTCCACCACCCCTTCCGG CAAACTCCAGTCAAAACCCTGCCCANCACTCCCAAAAAGTACTGCACCATAACCATCCCCATC AACTTCCCCGCGCAAATCNTCCTCCCGGNGCCGAACGGTATGANTNAANATCGTNNCCCANNG GANNAATCCTCGTCNNCTTNCNTTAGATANANNNNCTCNGNACNGTACTANTNNGGTAACNCTGA CTCANATGGGTNTCTCNNGATGTTTNNNNNNNTTTNNANCTAAAGNNCTTGTCGGNNNCGCGG AANTCNNGCCNGTCTNANTCCACCTTNTATAACTCTANGAGAGTTATNANNNGANCTCACGCA AATNAGTTCTNCTTNANGNCCCAAACNNTNNATGNANNNNACGACTTTACCTACCNGTATNAN ANCTTTTAACATNNCGNNCCANCCTGNATGCTCNTCTTNANNCNANANANNACTNNGAACTGT GTTNTANNATTATNTNNTANTNTTCTNTNNNNNNNCANNTTANCATTATNTNCTNNNNNNNNNTT NAAANCAANNATTNTNNNNANCATTCTNNTAATTAANATAATNCTCCNCGCNNNNNTTTCNAGANN NTNNCCNNNTNNNANNATTTCTTANATTTTTTANNNAATANNNNNTTTTTTNTNNNN



# BLAST Results


- Forward Read
  - 647 bp
  - 100% match to target sequence
- Reverse Read
  - 241 bp
  - 92% match to target sequence





# pSB1A3 & pSB1A7

- High copy number plasmids carrying ampicillin resistance
- pSB1A3 is 2,157 bp
- pSB1A7 is 2,431 bp





# Verifying Plasmid Identity

- Isolate plasmid DNA
  - Glycerol stocks of *E. coli* DH5alpha
  - Miniprep/GET buffer protocol from Openwetware.org
- Restriction digest using fast digest EcoRI
- Electrophoresis
  - do not give the accurate result which should be 2,157 bp for pSB1A3 and 2,431 bp for pSB1A7
  - Lots of undigested DNA
- Repeated the experiment and same result





Fig. 1. Electrophoresis of pSB1A3 & pSB1A7 plasmid DNA digested with fast digest EcoRI

# Other Extraction Attempts

- Extraction of plasmids from Parts Registry
  - We use registry kit of 2007 and paper punches from parts registry, spring 2008
  - Transform into NEB competent cell
  - Only the growth on control plates where we use pBLUE as a control
  - Parts not compatible with *E.coli* DH5alpha
  - Plasmids contain death gene called ccdB which kill all the competent NEB cells and occur no transformation.
- Extraction from different strain of *E. coli*, DB3.1
  - Digested with fast digest EcoRI, once again saw bands of incorrect size
  - But the control (pBLUE) works
  - Problem may be in the restriction enzyme we used, chemicals we prepared or in the plasmid itself
  - Repeat the same experiment, but different enzymes: fastdigest SpeI and XbaI and non fast digest EcoRI
  - Similar result
- Extraction from MIT samples
  - Digested with fast digest SpeI, fast digest XbaI, & non fast-digest EcoRI
  - pSB1A3 samples showed bands between 2,000 and 2,500 b.p
  - No clear bands appeared on the pSB1A7 samples



# Procedure Change

- Used different plasmid extraction procedure, GeneJet Plasmid Miniprep Kit to extract plasmid DNA from MIT samples
  - used SpeI and XbaI for digestion
  - Clearer bands, much less undigested DNA
  - Band sizes were opposite
    - pSB1A3 showed band size of approximately 2,500 bp
    - pSB1A7 showed band size of approximately 2,100 bp




Fig. 2. Electrophoresis of pSB1A3 & pSB1A7 plasmid DNA digested with SpeI & XbaI after Gene Jet extraction procedure

# Plasmid Verification using PCR

- PCR of pSB1A3 and pSB1A7 using VF2 and VR primers to amplify a segment of DNA
- VF2-VR= 316 bp in pSB1A3 and 590 bp in pSB1A7
- We observe VF2-VR is around 350 bp in pSB1A3 and around 900 bp in pSB1A7
- pSB1A7 is not target plasmid
- Also prove that the samples have been mislabeled while working in the laboratory only




Fig. 3. Electrophoresis of amplified fragment of pSB1A3 & pSB1A7 after PCR using VF2 & VR primers at a temp range of 48C – 52C

# Final Verification

- Extracted plasmid DNA using the Gene jet plasmid method
- Used SpeI and XbaI for digestion
- pSB1A3 showed bands at approximately 2,100 bp
- pSB1A7 showed bands at approx. 1,700 bp
- Concluded MIT samples of pSB1A7 contained “unknown” plasmid
- Also, MIT samples of pSB1A3 contained pSB1A3 plasmid described in registry





Fig. 3. Electrophoresis of pSB1A7 & pSB1A3 plasmid DNA samples digested with SpeI & XbaI

# Selection of promoters


- We selected
  - BBa\_R0040
  - BBa\_R0051
  - BBa\_R0010
  - BBa\_R0011
  - BBa\_I14032






# IDENTITY OF PROMOTERS


Bba\_R0040 Bba\_R0011 Bba\_R0051




# Double Digestion of pSB1A2 with Promoter with Spe1 & Pst1




# Double Digestion of PCR Product by Xba1 & Pst1



# Double digestion of cDNA with Xba1 & Pst1




# Double digestion of all RT PCR Samples with Xba1 & Pst1



# Problems

- Missed a restriction site
  - Fix: Reanalyze mRNA sequence with more than one program
- PCR amplification
- Maintaining Stocks at Correct Temperatures
  - Fix: Make sure everyone reads the label






# Special Thanks



Dr. Axel  
Schwekendiek



# Our Financier



=





# Remember

When all else fails try...

***RING OF POWER!***



How long your Prof.  
thinks it should take  
to do something



"Trivial"

=

How long it'll  
actually take you  
to do it



There goes your week.

"Easy enough"

=

Pull your hair out for  
a month.

"About a week"

=

Actually, this is pretty  
easy. He/she doesn't  
know there's technology  
that will do this for you  
now. Take the week off!

"Should keep you  
occupied for the rest  
of the term"

=

He/she will forget they  
asked you to do this by  
the end of the term.  
Don't even bother.

"This might make a  
good thesis topic"

=

Say hello to your  
thesis topic.

"Hmmm..."

=

Uh oh.