

# APPLICATION OF MICROTITER-BASED OXYGEN SENSORS TO ASSESS MICROBIAL COMMUNITIES

Sybil Sharvelle, Colorado State University


Jay Garland, Eric McLamore, M.K. Banks

# Overview



- Project description
- Description of oxygen sensor
- Applications
  - ▣ Surfactant biodegradation pathway
  - ▣ Spatial variation of functional groups of microorganisms in a biotrickling filter

# Biotrickling Filter for Graywater and Waste Gas Treatment



# Biological Conversions

- Surfactant removal from liquid phase (Pert Plus Kids)



- $NH_3$  gas solubilized to liquid phase where nitrification will occur



- $H_2S$  oxidation to liquid  $SO_4$




# Graywater Simulant

<i>Constituent</i>	<i>Concentration (mg/L)</i>	<i>NASA Potable Water Standard (mg/L)</i>
<b><u>Surfactants:</u></b>		
sodium laureth sulfate	73	0.5 <sup>a</sup>
disodium cocoamphodiacetate	62	0.5 <sup>a</sup>
polyalcohol ethoxylate	69	0.5 <sup>a</sup>
<b><u>Inorganics:</u></b>		
sodium	74	-
chloride	55	200
potassium	51	-
bicarbonate	9	-
ammonium	5	0.5
<b><u>Organics:</u></b>		
urea	87	0.5
lactic acid	58	0.5
acetic acid	38	0.5
formic acid	22	0.5
propionic acid	20	0.5


5

<sup>a</sup> – EPA Standard for foaming constituents

# Six Replicate Bench Scale Reactors




# Microtiter Based Oxygen Sensors



# Surfactant Biodegradation Pathways

PAE (polyalcohol ethoxylate)





(Steber and Wierich, 1985)

# Preferential Uptake of Metabolites

---


- Batch Assays


- Wastewater microbial consortium supplied with PAE as sole carbon source for at least 3 weeks
- Control population cultured in tryptic soy broth

- Previously identified metabolites of PAE biodegradation supplied to each culture in microtiter plate wells

- More rapid uptake or a higher response indicates previous exposure to that metabolite
  - Identify likely metabolites


# Response to PAE






(Steber and Wierich, 1985)

# Response of PAE Acclimated Micororganisms to Various Substrates





EG – ethylene glycol  
DEG – diethylene glycol  
TEG – triethylene glycol  
HEG – hexaethylene glycol  
OEG – octaethylene glycol



GOA – glyoxalic acid  
OA – oxalic acid  
GCA – glycolic acid  
FA – formic acid

# PEG Biodegradation



DEG and EG may accumulate

# Summary of PAE Biodegradation Assays

---

- Dominant pathway in studied systems is hydrophobe-phile scission as indicated by the formation of PEGs
- EG and DEG may accumulate
- Could be applied to bioreactor samples
  - ▣ Preliminary screening without intensive analytical methods development

# Spatial Variation of Functional Groups in Biotrickling Filters




# Experiment Description

---

- Coexistence of diverse microbial community (heterotrophs and autotrophs)
  - ▣ Biodegradation of surfactants
  - ▣  $\text{NH}_3 \rightarrow \text{NO}_3$
  - ▣  $\text{H}_2\text{S} \rightarrow \text{SO}_4$
- Biological samples scraped with scalpel from biotrickling filter packing material
  - ▣ Placed in microtiter plate wells
- Surfactants, ammonia and sulfide supplied to sample wells




# Reactor Sampling

BTF





- 9 samples collected in each replicate reactor
- 3 replicate reactors

# Gaussian Distributions



# Gaussian Distributions



# Summary of Spatial Variation Experiments

---

- Ammonia and sulfide oxidizing bacteria capable of coexisting with surfactant degraders in biotrickling filter
- DSCADA degraders not as well established as other surfactant degraders
- Ammonia and sulfide oxidizing bacteria located predominately at top of BTFs where contaminated gas entered

# Overall Summary

---

- Oxygen sensors are a useful tool to elucidate information about biological processes in engineered and natural systems
  - ▣ Rapid analysis of a large number of samples

# Questions